

DreamStation

Auto **BiPAP**

Bi-level therapy for obstructive sleep apnea patients: a comfortable alternative to CPAP

DreamStation Auto BiPAP features

- · Auto BiPAP, fixed BiPAP and fixed CPAP therapy modes
- · Adjustable pressure support (PS) : PS min-PS max
- Bi-Flex Comfort, C-Flex and C-Flex+
- · Advanced respiratory event detection*
- Smart Ramp : up to 45 min
- Mask fit check
- On-board memory
- Detailed compliance report
- Oximetry module (optional) to monitor therapy efficacy

Humidification

- Active humidification controlled by the room environment
- Pre-heat humidifier capable
- Heated tube (optional)

DreamStation Auto BiPAP benefits

- Daily Progress Feedback displays a simple trend of patients' nightly use to help track progress on the path to adherence.
- **EZ-Start** : automatic, personalized adjustments to CPAP pressure can help patients gradually acclimate to their prescribed level of therapy.
- **Connectivity** and remote monitoring options including a standard integrated Bluetooth module.
- **Performance Check** simplifies in-home device evaluation for easy remote troubleshooting.
- Patient' data management through detailed compliance and waveform reports.

*BiPAP auto detects and logs the following events: Open airway apnea (CA), Obstructed airway apnea (OA), Hypopnea (H), respiratory events related arousal (RERA), flow limitation (FL), vibratory snore (VS), large leak (LL), periodic breathing (PB).

DreamStation Auto BiPAP specifications

Part number	Description
INX700S15	DreamStation Auto BiPAP with SD card
INX700H15	DreamStation Auto BiPAP with humidifier and SD card

Specifications

Therapy modes	Auto BiPAP, fixed BiPAP, fixed CPAP
Therapy features	EZ-Start
Flex pressure relief	0 to 3
Pressure range	4-25 cm H ₂ O
Smart Ramp time	0 to 45 min (5-min increments)
Starting ramp pressure	Patient adjustable
Humidification	Heated humidification: fixed, adaptive, heated tube modes Heated tube (optional)
Device controls	Color LCD/control wheel/push button
LCD backlight control	Ambient light sensor
Filters	Reusable pollen and disposable ultra-fine
Data storage capacity (mininum)	Display: compliance 1 day, 7-day, 30-day averages SD card > 1 year On-board: up to 6 months
Altitude compensation	Automatic
Warranty	2 years
Dimensions	15,7 (L) x 19,3 (W) x 8,4 cm (H) (without humidifier) 29,7 (L) x 19,3 (W) x 8,4 cm (H) (with humidifier)
Weight	1,33 kg (without humidifier), 1,98 kg (with humidifier) but including power supply
Power supply	80W/12V/6,67A
Electrical requirements	100-240 VAC, 50/60 Hz, 2.0-1.0 A

Sleep therapy accessories	
Humidification	Part number
DreamStation heated humidifier	EUXH
Water chamber	1122520
Tubing	Part number
Heated tube (15 mm)	HT15
Standard tube (15 mm)	PR15
Filters	Part number
Pollen (foam) filter, reusable (1/pack)	1122446
Ultra-fine filter, disposable (1/pack)	1122447
Ultra-fine filter, disposable (2/pack)	1122518
Ultra-fine filter, disposable (6/pack)	1122519
Accessories	Part number
DreamStation carrying case	1121162
DreamStation link module/SpO ₂	1120293
DreamStation SpO ₂ cable	1121694
Reuseable finger sensor (clip)	936
Oximetry finger sensor, Adult, Flex	953
Flexiwrap, sens tape, Adult (25/pack)	954A
Reuseable finger sensor (soft), medium	1062050
Reuseable finger sensor (soft), small	1062051
DreamStation, 80W power supply	1118499
Connectivity	Part number
DreamStation cellular/SpO2 modem	100621C



Philips Respironics 1010 Murry Ridge Lane Murrysville, PA 15668 USA

DreamStation Wifi/SpO, modem

100721W

www.philips.com/respironics