

# Reusable Soft Sensors

Now with  
2-year warranty



## The Universal Sensor

For spot checks or continuous monitoring, Nonin's soft sensors perform like two sensors in one, with patient comfort in mind.

- Quick and easy spot check application
- Comfortable fit for continuous monitoring
- Durable design delivers consistent performance
- Suitable in many medical settings (EMS, Hospital, Sleep)
- Easy-to-clean with isopropyl alcohol
- Latex-free
- Two-year warranty

## Sizes Adult – Infant



**8000SL - Large**  
Digit thickness 12.5 - 25.5 mm.  
(.5 and 1 inch)



**8000SM - Medium**  
Digit thickness 10 - 19 mm.  
(.4 and .75 inches)



**8000SS - Small**  
Digit thickness 7.5 - 12.5 mm.  
(.3 and .5 inches)



## PureLight® Accuracy

PureLight® technology means accurate readings are assured.

**Blood Oxygen Saturation**  
(%SpO<sub>2</sub> A<sub>rms</sub>\*) 70 - 100% ±2 digits for adults

**Pulse Rate**  
±3% ±1 digit

[www.nonin.com](http://www.nonin.com)

800.356.8874 +1.763.553.9968

Please contact your authorized distributor for more information.

\* A<sub>rms</sub> encompasses 68% of the population.

