

**AeroChamber<sup>®</sup>**  
Plus\*


Valved Holding Chamber



## Proven Performance

Getting the medication to your patients' lungs is priority one - the **AeroChamber Plus\*** Valved Holding Chamber (VHC) is designed to do just that.

- Proven dosing performance across a wide range of formulations
- The only full line of aerosol delivery devices designed specifically for each patient
- **EZ-Flow\*** Exhalation Valve offers lower resistance to exhaled flow and directs the flow away from the patient to prevent irritation



*ComfortSeal<sup>\*</sup> Mask provides a comfortable, secure fit*



Trudell Medical International™

[www.trudellmed.com](http://www.trudellmed.com)

**EASY TO USE... EASY TO LEARN... MOST PREFERRED.<sup>1</sup>**



**AeroChamber<sup>®</sup> Plus\***  
Valved Holding Chamber



### ComfortSeal<sup>\*</sup> Mask

Latex-free, soft, comfortable, easy to clean and provides a secure seal across a wide age range



### EZ Flow<sup>\*</sup> Exhalation Valve

Robust design allowing a low resistance to exhaled flow



### Flow Dynamic Design

Increases suspension time for particles, improves respirable dose of medication

**EASY TO USE... EASY TO LEARN... MOST PREFERRED.<sup>1</sup>**

<sup>1</sup> NOP Healthcare, UK, March, 2000.



Trudell Medical International<sup>®</sup>

Trudell Medical International  
725 Third Street, London,  
Ontario, Canada, N5V 5G4  
Phone: +1-519-455-7060  
Fax: +1-519-455-7858  
[www.trudellmed.com](http://www.trudellmed.com)

Authorised European Representative:  
Trudell Medical International Europe Ltd.  
BioCity Nottingham, Pennyfoot Street  
Nottingham, NG1 1GF, UK  
Tel: +44 115 912 4380  
Fax: +44 115 912 4289

